DENSORAPID D

Handheld Densitometer DENSORAPID D



Measuring optical density of X-Ray films with high accuracy

DENSORAPID D handheld Densitometer

For precise determination of optical density of X-Ray films, DENSORAPID D is equipped with an absolute innovative digital measuring technique. The illuminated display allows easy reading even in darkened rooms and the calibration process is highly simplified by the sophisticated software.

Fault-free measurements even in slightly inclined positions are guaranteed by the ergonomic hand probe with a soft tip avoiding scratches on the film.

The battery is charged with a standard USB charger and ensures an operating time of up to 30 hours.

Specification:

Measuring range: Resolution: Accuracy: Repeatability: Operating time: USB charger: Dimensions: Weight: 0.00 ... > 5.00 D 0.01 D +/- 0.02 D +/- 0.01 D 30 hours with fully charged battery 100...240V Micro-USB 135 x 73 x 29 mm 0.25 kg

Scope of supply:

DENSORAPID D in transport case complete with magnetic measuring mask 10x24 cm / Ø 10 mm, USB charger, KOWOSTEP D density step tablet, calibration certificates Article no. 13 40540

Accessories: USB charger Article No.: 13 40546 Transport case Article No.: 13 40547

Measuring mask, magnetic, 10x24 cm / Ø 10 mm *Article No.: 13 40104* Measuring mask, magnetic, 10x48 cm / Ø 10 mm *Article No.: 13 40084*

KOWOSTEP D density step tablet Article No.: 13 40207 Recalibration KOWOSTEP D, 3 years interval Article No.: 13 40208 BAM density step tablet Article No.: 13 40205 Recalibration BAM density step tablet, 3 years interval Article No.: 13 40206

Densitometer mounting device for KOWOLUX X film viewer, universal Article No.: 13 40087

















Solinger Strasse 186 D-40764 Langenfeld Fon: +49 2173 - 223 83 Fax: +49 2173 - 223 35

3 www.kowotest.de 5 info@kowotest.de